

Figure 1

## **System Architecture**









Figure 4



49, THE ENTIRE SKIN SURFACE








Figure 4, Cont.









4e. SHOWLDERS, CHEST, BACK & ARMS







Figure 4, Cont.




4g. CHEST, BACK, TRRSO, \$ LEGS



## Control Alg rithm - Co ling



## Temperature Gradient [Cooling]



 $\Delta T$  Temperature Gradient  $= |T_{Core} \cdot T_{Interface}|$ 

Heat Transfer is the Driving Force in: at the


Thermal Interface

• Cooling:  $T_{Interface} < T_{Core}$ 

• Warming: T<sub>Interface</sub> > T<sub>Core</sub>

Figure 7

#### TInterface affects Vasoconstriction & Vas dilati n (as measured by Blood Flow)



For each individual,

• Vasoconstriction [VC] occurs below a certain Temp range

· Vasodilation [VD] occurs

above that Temp range


Blood Flow can be measure by:

- · Laser Doppler
- Bio-Impedance
- Light Absorption (Pulse Oximetry)

Figure 8

## Heat Transfer = $f(Temp Grad \times Blood Flow)$


Figure shows Temp Grad & Blood Flow vs.  $T_{\text{interface}}$  superimposed



# Maximum Heat Transfer occurs @

The lowest T<sub>Interface</sub> where Vasodilation occurs

#### Hysterysis:



T<sub>Interface</sub> [°C]


The transition between Vasoconstriction and Vasodilation is NOT Identically Reversible...

The transition occurs at a different temperature range depending on the initial condition

Typically, the transition from:

VC → VD occurs at a T<sub>Interface</sub> range above  $VD \longrightarrow VC$ 


### If Vas dilati n is initially det cted



- A Blood Flow Sensor detects VD, T<sub>interface</sub> = T<sub>set</sub>
- B System controller decreases T<sub>interface</sub> until VC detected
- C T<sub>interface</sub> increases above transition temp range, VD occurs
- D System controller decreases T<sub>interface</sub> to T<sub>Max Heat Transfer</sub>

T<sub>Max Heat Transfer</sub> < T<sub>set</sub>

## If Vasoc nstricti n is initially detect d



- A Blood Flow Sensor detects VC, T<sub>Interface</sub> = T<sub>set</sub>
- B System controller increases T<sub>Interface</sub>
- © T<sub>interface</sub> increases above transition temp range, VD occurs
- D System controller decreases T<sub>Interface</sub> to T<sub>Max Heat Transfer</sub>

T<sub>Max Heat Transfer</sub> > T<sub>set</sub>